

10/501262  
DT04 Rec'd PCT/PTO 09 JUL 2004

IN THE UNITED STATES RECEIVING OFFICE (RO/US)  
Designated/Elected Office (DO/EO/US)

U.S. National Stage of

International Application No.: PCT/GB03/00064

International Filing Date: 10 January 2003

Priority Date Claimed: 11 January 2002

Applicant: William Malcolm Charles Rosenberg

Title: Methods of Detecting HCV Genotype 1 (HCV-1) by  
Using Primers Specific for the 5' Non-Coding Region  
(NCR) of the HCV Genome

Attorney's Docket No.: 3769.1000-000

Date: 09 July 2004

EXPRESS MAIL LABEL NO. EV215891406US

11/PRZS

10/501262

DT04 Rec'd PCT/PTO 09 JUL 2004

PCT/GB03/00064

WO 03/057915

1

## METHODS

This invention relates to methods for detecting viral genotypes, and in particular Hepatitis C virus (HCV) genotypes. The invention therefore finds application in  
5 determining prognosis, and in the selection of treatment regimes, for patients infected with HCV.

The hepatitis C virus is a positively stranded RNA virus that exists in six genetically distinct genotypes. These are designated Type 1, 2, 3, 4, 5 and 6, and their full length  
10 genomes have been reported (Genbank/EMBL accession numbers Type 1a: M62321 (Choo), AF009606, AF011753, Type 1b: AF054250, D13558, L38318, U45476, D85516; Type 2b: D10988; Type 2c: D50409; Type 3a: AF046866; Type 3b: D49374; Type 4: WC-G6, WC-G11, WG29 (Li-Zhe Xu *et al*, *J. Gen. Virol.* 1994, 75:2393-98), EG-21, EG-29, EG-33 (Simmonds *et al*, *J. Gen. Virol.* 1994, 74: 661-  
15 668)). Viruses in each genotype exist as differing "quasispecies" that exhibit minor genetic differences. In the UK and USA, the vast majority of infected individuals are infected with genotype 1, 2 or 3 HCV. In the UK, the only other genotype encountered is genotype 4 which is identified at low prevalence in some immigrant populations. HCV infection affects approximately 0.5% of the population in the  
20 United Kingdom, 1.8% of the population in the USA and 3% of the population of the world. In over 85% of infected people, HCV causes a lifelong infection characterised by chronic hepatitis that varies in severity between individuals.

The only effective therapy that is currently available is interferon alpha (IFN- $\alpha$ ) used  
25 alone or in conjunction with ribavirin ("combination therapy"). Randomised controlled trials have shown that combination therapy is highly effective in the treatment of HCV infection in patients previously treated with interferon alone and in patients never previously treated with interferon (Davis *et al*, *NEJM*, (1998), 339(21): 1493-99; Poynard *et al*, *Lancet* (1998) 352(9138): 1426-32).

Subgroup analysis has revealed that the major determinant of prognosis and response to treatment is HCV genotype (Poynard *et al*, *Lancet*, 1998, 352(9138): 1426-32). For patients with genotype 2 or 3 infections treated for six months, response rates to treatment for six months are equivalent to those seen in patients with genotype 1 infections for one year. There were insufficient data in the trials on response rates for patients with genotype 4, 5 and 6 infections to permit investigators to determine response rates in patients with these infections. In mixed genotype infections, genotype 1 predominates.

Combination therapy is expensive (in the UK, it costs approximately £8,000 for one year's treatment). Thus, the identification of patients with non-genotype 1 infections allows treatment to be stopped after 6 months with a consequent cost saving, and the avoidance of unnecessary adverse effects of the drugs associated with further treatment. In addition, the most serious outcome of HCV infection is hepatocellular carcinoma. It is therefore desirable to be able to identify with which genotype(s) of HCV a patient is infected.

Known HCV genotyping methods rely on agarose gel electrophoresis of restriction endonuclease digested polymerase chain reaction amplification products (RFLP) or antibody detection of anti-HCV antibodies. Enzymatic digestion can take a relatively long time to reach completion making the RFLP method time consuming. It also involves several different stages, making it unamenable to automation or rapid throughput. In addition, inaccurate genotyping can occur owing to (i) the fact that restriction endonucleases do not cut with 100% efficiency, and (ii) single base pair mutations in quasispecies may result in failure of digestion of the PCR products. Antibody based methods are costly and relatively insensitive. Most significantly, they are incapable of differentiating between patients who have cleared HCV infection and those with persisting infection. They thus require the performance of a virological test in parallel to confirm current infection.

US Patent No. 5851759 discloses methods of genotyping HCV in which HCV RNA is isolated and cDNA is synthesised from this RNA. The cDNA is then subjected to PCR using primers which flank the E1 region of the HCV genome. The products of PCR are analysed using the Heteroduplex Tracing Assay (HTA). The E1 region was selected for analysis because it is considered to be the most heterogeneous region of the HCV genome. Thus, amplification of this region is likely to give a different product for each HCV strain in a sample. These products can then be analysed to indicate the strains that are present in the sample.

10

In view of the different effect of combination therapy on infection with HCV 1 as compared to infection with other HCV subtypes, there exists a need to be able to detect HCV genotype 1.

15

According to the present invention, there is provided a method of detecting HCV genotype 1 (HCV-1) in a sample, comprising:

subjecting the sample to an amplification reaction using at least one primer which anneals specifically to the 5' noncoding region (5' NCR) of the HCV-1 genome; and:

20

detecting the product of the amplification reaction.

The HCV genome consists of 5' and 3' noncoding (or untranslated) regions (NCR) that flank a single long open reading frame (ORF). This ORF encodes for three structural proteins at the amino-terminal end and for six non-structural (NS) proteins at the carboxyl-terminal end. The structural proteins are represented from the nucleocapsid (core; C) proteins and two glycoproteins, envelope 1 (E1) and envelope 2 (E2). The non-structural proteins are named NS2, P7, NS3, NS4a, NS4b, NS5a, NS5b. The 5' NCR is the most highly conserved part of the HCV genome, whereas the sequence of the two envelope proteins (E1 and E2) is highly variable among different HCV isolates. The present invention lies in the surprising finding that parts of the 5' NCR

30

are conserved between HCV-1 quasi species but not between other genotypes. A preferred region is the region between residues -134 and -118 of the 5'NCR of HCV-1. Accordingly, primers for amplification reactions can be designed which anneal specifically to those parts of the 5' NCR that are conserved between HCV-1 quasi species but not between other genotypes.

The amplification reaction may be PCR (see for example U.S. Patents Nos. 4,683,195 and 4,683,202, and Innis *et al*, editors, PCR Protocols, (Academic Press, New York, 1989; Sambrook *et al*, Molecular Cloning, Second Edition, (Cold Spring Harbour Laboratory, New York 1989)). PCR will generally be used when HCV RNA has been isolated and converted, preferably by reverse transcription, to cDNA. Preferably, PCR is carried out using Taq DNA polymerase, e.g. Amplitaq™ (Perkin-Elmer, Norwalk, Conn.). Taq polymerase can also be obtained from MBI Fermentas, Perkin Elmer, Boehringer Mannheim and Beckman Instruments. An equivalent, preferably thermostable, DNA polymerase may also be used in the method of the present invention, such as Tfl (*Thermus flavus*) polymerase (Gut *et al*, *Virol. Methods* 77(1): 37-46 (1999)).

Alternatively, the amplification reaction may be RT-PCR (Yajima *et al*, *Clin. Chem*, 44(12): 2441-2445 (1998); Martell *et al*, *J. Clin. Microbiol.*, 37(2): 327-332 (1999); Gut *et al*, *Virol. Methods* 77(1): 37-46 (1999); Predhomme *et al*, *Leukemia*, 13(6): 957-964 (1999)), in which the HCV RNA is reverse transcribed into cDNA which is then subjected to PCR amplification.

As is well-known, PCR involves the extraction and denaturation (preferably by heat) of a sample of DNA (or RNA). A molar excess of oligonucleotide primers is added, along with a polymerase, which may be heat-stable, and dNTPs for forming the amplified sequence. The oligonucleotide primers are designed to hybridise to opposite ends of the sequence which it is desired be amplified. In the first amplification round, the polymerase replicates the DNA to produce two "long products," which begin with

the respective primers. The total DNA, which includes the two long products and the two original strands, is then denatured and a second round of polymerisation is carried out (for example, by lowering the temperature). The result of the second round is the two original strands, the two long products from the first round, two new long  
5 products (produced from the original strands), and two "short products" produced from the long products. These short products have the sequence of the target sequence (sense or antisense) with a primer at each end. For each additional amplification round, the number of short products grows exponentially, each round producing two additional long products and a number of short products equal to the sum of the long  
10 and short products remaining at the end of the previous round.

Oligonucleotide primers can be synthesised by a number of approaches, e.g. Ozaki *et al*, *Nuc. Acids Res.* **20**: 5205-5214 (1992); Agrawal *et al*, *Nuc. Acids Res.* **18**: 5419-5423 (1990) or the like. Conveniently, the oligonucleotide probes are synthesised on  
15 an automated DNA synthesiser, e.g. an Applied Biosystems, Inc, Foster City, California model 392 or 394 DNA/RNA synthesiser using standard chemistries such as phosphoramidite chemistry (Beaucage and Iyer, *Tetrahedron* **48**: 2223-2311 (1992), US Patent Nos. 4980460, 4725677, 4415732, 4458066 and 4973679). Alternative chemistries, including non-natural backbone groups such as phosphorothioate and  
20 phosphoramidate, may also be employed, provided that the hybridisation efficiencies of the resulting oligonucleotides are not adversely affected. The precise length and sequence of the DNA primers will depend on the target polynucleotide to be amplified. Preferably, the length of the DNA primers is in the range 10 to 60 nucleotides and more preferably in the range 15 to 30 or 25 nucleotides.

25 In a preferred embodiment, the method of the present invention further comprises a step in which HCV infection generally is detected. In this step, the sample is subjected to an amplification reaction using primers which anneal to a region of the 5'NCR which is conserved between all HCV genotypes and the product of this  
30 amplification reaction is detected. Preferably, the amplified region spans the region

which is conserved in HCV-1. This additional step allows verification that HCV is present in the sample and so treatment can be tailored accordingly if no HCV-1 is detected.

- 5 Detection of the products of amplification may be carried out using the well-known technique of agarose gel electrophoresis. If HCV-1 is present in the sample, the amplification reaction produces a product (of known size) which can be detected on the agarose gel. An HCV-1 specific primer suitable for this is the DNA oligonucleotide:

10

5' CCI CTC AAT GCC TGG AG 3' (*Spec-1*)

*Spec-1* may be used together with any universal HCV primer. One example is:

15

Reverse: 5' GCA GTA CCA CAA GGC CTT TCG C 3' (*UTR-R2*)

The invention provides such primers alone and in combination.

- 20 If it is desired to carry out an amplification reaction additionally to detect whether HCV in general is present in the sample, a pair of primers universal for all HCV genotypes can be used. This will produce a product of a known size which is different to that for the HCV-1-specific amplification. If HCV-1 and other HCV genotypes are present, the gel will show both products. Suitable primers for this are:

25

Forward: 5' CGT CTA GCC ATG GCG TTA G 3' (*UTR-L2*)

Reverse: (*UTR-R2*)

30

Additionally, a preliminary amplification may be carried out on the sample to isolate HCV material. This preliminary amplification can be carried out using two primers universal for all HCV genotypes. Suitable primers for this are:

Forward: 5' GGA ACT ACT GTC TTC ACG C 3' (*UTR-L1*)

Reverse: 5' ACG GTC TAC GAG ACC TC 3' (*UTR-R1*)

5 Detection of the products of amplification may alternatively be carried out using  
fluorescent analysis. If genotype 1 HCV is present, amplification occurs and this is  
detected by hybridisation of a quenched fluorescently labelled probe. With successive  
rounds of amplification, an increasing amount of the fluorescent probe becomes  
incorporated into PCR products and the quencher is liberated resulting in an increasing  
10 amount of fluorescence being detectable on the reaction vessels. Such fluorescence  
detection makes this approach ideally suited to automation and high throughput  
processes.

One well-known fluorescent technique is the so-called "TaqMan" method (Holland, *et*  
15 *al*, *Proc. Nat. Acad. Sci*, **88**: 7276-80, (1985)). This method uses an oligonucleotide  
probe complementary to a segment of the target polynucleotide, i.e. between the  
sequence defined by the PCR primers, which probe is labelled with two fluorescent  
molecules or fluorophores. The emission spectrum of one of the molecules overlaps  
with the excitation spectrum of the other and, as a result, the emission of the first  
20 fluorophore (the "reporter") is quenched by the second (the "quencher") when they are  
in sufficient proximity.

The oligonucleotide probe is hybridised to the target polynucleotide downstream from  
a primer for a polymerase having a 5'-3' exonuclease activity. As the polymerase  
25 extends the primer as part of the PCR reaction, the oligonucleotide is digested,  
releasing one of the "quencher" and "reporter" molecules and causing the distance  
between these molecules to become such that the fluorescent emission is no longer  
quenched by energy transfer. Thus, a fluorescent signal is generated, providing real-  
time monitoring of amplification.



The oligonucleotide probe may be made in the same way as the oligonucleotide primers mentioned above and may have the same variations in backbone and so on, providing that hybridisation to the target polynucleotide is not compromised and/or cleavage efficiency of the exonuclease are not adversely affected. Preferably, the length of the oligonucleotide probe is in the range of 10 to 60 nucleotides, and more preferably 18 to 30 nucleotides. The precise length and sequence of the oligonucleotide probe depends at least in part on the nature of the target polynucleotide to which it binds and will be varied to achieve appropriate annealing and melting properties for a particular target polynucleotide. In addition, the binding location of the probe may be varied according to the nature of the target polynucleotide.

Preferred primers for this reaction are *UTR-R2* and *Spec1* mentioned above. A preferred probe has the following sequence and is specific to HCV-1:

5' FCG CIA CCC AAC ICT ACT IGG CTA GT 3' (*L1*)

where F=6-FAM, 3'-T+TAMRA

Alternatively, the technique such as the TaqMan technique can be used to detect HCV-1 using one or more universal primers and an HCV-1 specific probe. Thus, in a further aspect, the present invention provides a method of detecting HCV genotype 1 (HCV-1) in a sample, comprising:

subjecting the sample to an amplification reaction using at least one primer which anneals to the genome of HCV, a polymerase having a 5'-3' exonuclease activity and an oligonucleotide probe, which probe anneals specifically to the 5' noncoding region (5' NCR) of the HCV-1 genome and which incorporates a modified nucleotide having a fluorescent characteristic which is modified by one or more neighbouring nucleotides; and

detecting a change in fluorescence as the oligonucleotide probe is degraded by the exonuclease activity of the polymerase as the polymerase extends the primer and modification of the fluorescent characteristic of the modified nucleotide is reduced.

5 The modified nucleotide is preferably a "reporter" molecule and the neighbouring nucleotide(s) is/are preferably a "quencher" molecule. The primer(s) may be universal and is/are preferably *UTR-L1* and *UTR-R1* or *UTR-L2* and *UTR-R2*. Alternatively, the probe(s) may be capable of annealing specifically to the 5' NCR of the HCV-1 genome. The probe may be *L1* as defined above.

10

A further alternative method of detecting the amplification products is using one or more molecular beacon probes. These are PCR primers comprising a sequence attached to a hairpin loop sequence which contains both a fluorescent label and a quencher molecule. The hairpin brings the fluorescent molecule in opposition to the quencher so that the primer does not fluoresce under normal conditions. As PCR results in copying of the probe sequence, the hair-pin opens and the fluorescent beacon and quencher become separated. The resulting PCR product fluoresces.

15

In a preferred embodiment, the following PCR primers are used:

20

Forward: *MBP-LR-1*

5'FCA CCT TCA CCC TCA GAA GGM GCC GCT CAA TGC CTG GAG3'

(F=FAM; M=MeREDdU and U=Uracil)

25

Reverse: *UTR-R2*

*MBP-LR-1* is a forward HCV type 1 specific primer and, as mentioned above, *UTR-R2* is a universal reverse primer. Any universal HCV primer can be used together with *MBP-LR-1*.

30

If it is desired to detect the presence of other HCV genotypes in the sample, an additional amplification can be carried out using at least one molecular beacon probe which is universal for all HCV genotypes. A suitable probe in this regard is *MBP-LR-ALL* which is a forward primer and has the following sequence:

5

5'FCACCTTCACCCTCAGAAGGMGCGUCTAGCCATGGCGTTAG3'  
(F=FAM; M=MeREDdU and U=Uracil)

10 *MBP-LR-ALL* can be used together with any universal reverse primer; one suitable example is *UTR-R2*.

Fluorochromes that emit light at differing wavelengths can be incorporated into the different primers so that products amplified with different primers can be distinguished in the same reaction vessel. Alternative hairpin sequences on the beacon probes can easily be designed by the skilled person.

Alternatively, the amplification products can be detected using DNA hybridisation employing enzyme linked methods to confirm hybridisation, such as using horse radish peroxidase on sequence specific probes and an appropriate substrate.

20

According to a further aspect of the present invention, there is provided a kit for detecting HCV genotype 1 (HCV-1) in a sample, comprising: at least one primer which anneals specifically to the 5' noncoding region (5' NCR) of the HCV-1 genome.

25

The kit may further comprise the polymerase and an appropriate primer or set of primers. In addition, it may comprise additional reagents that are necessary for performing the invention, such as dNTP mixtures, buffers, molecular size standards, wax beads and the like. The reagents may be provided in pre-measured amounts so as to simplify the

performance of the method. Typically, the kit may also contain detailed instructions for carrying out the method.

5 It will be appreciated by those skilled in the art that the sequences of the primers and probes described above can be modified without affecting their activity, i.e. their ability to act as primers or probes in the methods of the present invention. Such modified probes and primers are included in the present invention, provided that they have substantial identity with the probes mentioned above.

10 When comparing nucleic acid sequences for the purposes of determining the degree of homology or identity, one can use programs such as BESTFIT and GAP (both from the Wisconsin Genetics Computer Group (GCG) software package). BESTFIT, for example, compares two sequences and produces an optimal alignment of the most similar segments. GAP enables sequences to be aligned along their whole length and finds the  
15 optimal alignment by inserting spaces in either sequence as appropriate. Suitably, in the context of the present invention when discussing identity of nucleic acid sequences, the comparison is made by alignment of the sequences along their whole length.

20 Preferably, sequences which have substantial identity have at least 75% sequence identity and more preferably at least 90 or at least 95% sequence identity with said sequences. In some cases, the sequence identity may be 99% or above.

25 Desirably, the term "substantial identity" indicates that said sequence has a greater degree of identity with any of the sequences described herein than with prior art nucleic acid sequences.

30 Preferred features of each aspect of the invention are as for each of the other aspects *mutatis mutandis*. The prior art documents mentioned herein are incorporated to the fullest extent permitted by law. The abbreviations for nucleotides used herein are standard unless otherwise indicated.

The present inventors have invented a polymerase chain reaction amplification (PCR) based test that permits the detection of all types of HCV infection and the specific identification of HCV genotype 1 infections.

5

The present invention will now be described further in the following examples. Reference is made to the accompanying drawings:

Figures 1a and 1b show 2% agarose gels of PCR products stained with ethidium bromide to show bands 88 base pairs in size (genotype 1 specific band) and 217 base pairs (all HCV genotypes).

10

Figures 2a-f show printouts of PCR runs using fluorescent probe detection and Taqman technology. Fluorescence intensity is represented on the y axis and the number of PCR cycles elapsed is shown on the x axis.

15

Figures 3a and 3b are printouts of PCR runs using fluorescent probe detection and Taqman technology. Fluorescence intensity is represented on the y axis and the number of PCR cycles elapsed is shown on the x-axis.

20

Figure 4 shows printouts of PCR runs using molecular beacon probe analysis and fluorescence detection. Fluorescence intensity is represented on the y axis and the number of PCR cycles elapsed is shown on the x axis.

25

Figure 5 is a diagram illustrating the positions of the various primers and probes used in certain embodiments of the present invention within the 5'NCR of the HCV-1 genome.

## Examples

### *Example 1*

- 5 From a 200µl sample of serum, HCV RNA was extracted using mRNA extraction kit (QIAGEN Ltd, Crawley, West Sussex) as per manufacturer's instructions. 11.5µl of RNA solution as extracted above was added to 0.5 µl of random primers (Promega UK Ltd, Southampton) in a 1.5 µl eppendorf test tube and incubated at 90 °C for 5 minutes in a heating block and then snap cooled on ice.
- 10 cDNA was then prepared by reverse transcription using the following reaction mix:
- 5 µl of 5x SuperScript buffer (Life Technologies (GIBCO BRL) Paisley, Scotland); 2.5 µl of 100 mM DTT; 2.5 µl of 1mg/ml Bovine Serum Albumin acetylated (10 mg/ml) (Promega, UK Ltd, Southampton); 1.25 µl of 10mM dNTP solution
- 15 (Amersham Pharmacia, Buckinghamshire); 1.25 µl Superscript Moloney Murine Leukaemia Virus reverse transcriptase (Life Technologies (GIBCO BRL) Paisley, Scotland); 0.5 µl of RNasin ribonuclease inhibitor (Promega, UK Ltd, Southampton).
- 20 13 µl of this solution was added to each tube containing RNA/random primers. The mix was incubated at 42 °C for 60 minutes. The reaction was stopped by heating to 95 °C for 5 minutes. The cDNA was then stored at -20 °C.
- 25 1µl of the resulting complimentary DNA (cDNA) was suspended in 20µl of a first round reaction mix. This mix contains 1.5mM MgCl<sub>2</sub>, 0.2mM deoxy nucleotide triphosphates (guanine, adenosine, threonine and cytosine) (Amersham Pharmacia, Amersham, UK), 0.1µl of Taq polymerase (Promega UK Ltd, Southampton) and 500nM each of primers *UTR-L1* and *UTR-R1* (from Cuachem Ltd, Glasgow).

*UTR-L1* is the forward primer and is complementary to residues -273 to -291 of the 5'UTR of HCV1. It has the sequence:

5'GGA ACT ACT GTC TTC ACG C3'

5

*UTR-R1* is the reverse primer and is complementary to residues -6 to -22 of the 5'UTR of HCV1. It has the sequence:

5'ACG GTC TAC GAG ACC TC3'

10

The PCR reaction conditions were as follows: 94°C for 1 minute for 1 cycle; 94°C for 30 seconds; 55°C for 20 seconds; 72°C for 20 seconds for 25 cycles.

1µl of this reaction was then added to 20µl of a second round reaction mix. This mix is the same as the first round reaction mix, except that primers *UTR-L1* and *UTR-L2* are replaced by 500mM of primers *Spec-1* and *UTR-R2*, and 100mM of *UTR-L2*.

15

*UTR-L2* is a forward primer which is universal for all HCV genotypes and is complementary to residues -248 to -266 of the 5'UTR of HCV1. It has the following sequence:

20

5'CGT CTA GCC ATG GCG TTA G3'

*Spec-1* is a forward primer which is HCV type 1 specific and is complementary to residues -118 to -134 of the 5'UTR of HCV1. It has the following sequence:

25

5'CCI CTC AAT GCC TGG AG3'

*UTR-R2* is a reverse primer which is universal for all HCV genotypes and is complementary to residues -50 to -71 of the 5'UTR for HCV1. It has the following sequence:

5                    5'GCA GTA CCA CAA GGC CTT TCG C3'

The reaction conditions were as follows: 94°C for 1 minute for 1 cycle; 94°C for 30 seconds; 58°C for 20 seconds; 72°C for 20 seconds for 25 cycles.

10        PCR products were then resolved by agarose gel electrophoresis (NuSieve 3:1, Flowgen Instruments, Lichfield, Staffordshire, UK) run for 40 minutes at 100 Volts in 2% agarose.

15        Samples of known HCV genotypes (determined by RFLP studies) were tested. The results are shown in Figures 1a and 1b. In Figure 1a, the lanes contain (from left to right):

1 – molecular weight markers  
2 – HCV negative control  
20        3 – HCV 3a  
          4 – HCV 2b  
          5 – HCV 1a  
          6 – HCV 1a  
          7 – negative control  
25        8 – HCV 1b  
          9 – HCV 1b  
         10 – water control  
         11 – HCV 1b  
         12 – HCV 2b



- 13 – negative control
  - 14 – HCV 1a
  - 15 – HCV 1b
  - 16 – HCV 3a
  - 5 17 – HCV negative control
  - 18 – 1a positive control
  - 19 – water control
  - 20 – molecular weight markers
- 10 In Figure 2a, the lanes contain (from left to right):
- 1 – molecular weight markers
  - 2 – negative control
  - 3 – HCV 1a control
  - 15 4 – HCV 3a control
  - 5 – HCV negative control
  - 6 – HCV 5
  - 7 – negative control
  - 8 – HCV 3a
  - 20 9 – negative control
  - 10 – HCV 2a
  - 11 – HCV 1b
  - 12 – HCV 1a
  - 13 – negative control
  - 25 14 – HCV 3a
  - 15 – HCV 1a
  - 16 – negative control
  - 17 – molecular weight markers

As can be seen, all HCV genotypes result in a 217 base pair band. Genotype 1 HCV generates an additional 88 base pair product.

*Example 2*

5

RNA was isolated and cDNA produced in the same manner as in Example 1.

Two separate 25µl reaction mixes containing 1µl of cDNA, 3.5mM MgCl<sub>2</sub>, 0.2mM each of dATP, dCTP, dGTP and 0.4mM dUTP, 0.25U of AmpliTaq Gold DNA Polymerase (Perkin Elmer), 300nM of *UTR-R2* and 300nM of *L-1* probe were prepared. One contained 300nM of *Spec-1* and the other contained 300nM of *UTR-L2*.

10 *L1* is a HCV genotype 1 specific probe (prepared by Scandinavian Gene Synthesis AB, Koping, Sweden) which is complementary to residues -68 to -93 of the 5'UTR of HCV1 and has the following sequence:

5'FCG CIA CCC AAC ICT ACT IGG CTA GT3'

20 where F=6-FAM,3'-T+TAMRA

The reaction conditions were as follows: 50°C for 2 minutes, 95°C for 10 minutes; for 1 cycle, 95°C for 15 seconds, 60°C for 60 seconds; for 60 cycles. Fluorescence was detected at 520nm.

25

Samples of known HCV genotypes (determined by RFLP studies) were tested. The results are shown in Figures 2a-f. The following table provides a key to the conditions:

	1	2	3	4	5	6
Primers	(UTR-L2 + UTR-R2)	(Spec-1 + UTR-R2)	(UTR-L2 + UTR-R2)	(Spec-1 + UTR-R2)	(UTR-L2 + UTR-R2)	(Spec-1 + UTR-R2)
A	-ve control	-ve control	sample	Sample	2a	2a
B	-ve control	-ve control	1a	1a	1b	1b
C	-ve control	-ve control	2b	2b	Sample	Sample
D	-ve control	-ve control	3a	3a	sample	Sample
E	-	-	1a	1a	3a	3a
F	1a	1a	5	5	1b	1b
G	sample	sample	1b	1b	3a	3a
H	1b	1b	1a	1a	2b	2b

Thus, it can be seen that HCV-1 can be differentiated from the other subtypes when the primer specific for HCV-1 is used (*Spec-1*) (see Figures 2b, d and f).

5 Similar results are seen in Figures 3a and 3b where the conditions are as follows:

	1 (UTR-L2 + UTR- R2)	2 (UTR-L2 + UTR- R2)	3 (UTR-L2)	4 (UTR-L2)	5 (SPEC-1 + UTR-R2)	6 (SPEC-1 + UTR- R2)	7 (SPEC-1 + UTR- R2)	8 (SPEC-1 + UTR- R2)
A	1a	3a	-ve control	-ve control	1a	3a	-ve control	-ve control

10 This approach relies on sequence specific hybridisation of a fluorescently labelled probe to genotype 1 HCV sequences. These sequences are generated by PCR using at least one primer specific for HCV-1. The probe has a quencher molecule linked to the fluorochrome. As PCR proceeds, the quencher molecule is removed by Taq polymerase from incorporated probes. Fluorescence of the resulting genotype 1 PCR

products can then be detected. The amount of fluorescence generated during the successive rounds of PCR correlates with the amount of PCR product.

### Example 3

5

RNA was isolated and cDNA produced in the same manner as in Example 1.

Two separate 25µl reaction mixes containing 1µl of cDNA, 3.5mM MgCl<sub>2</sub>, 0.2mM each of dATP, dCTP, dGTP and 0.4mM dUTP, 0.1 µl of Taq polymerase stoeffel  
10 fragment (Perkin Elmer), 0.1 µl of antiTaq antibody (Gibco BRL – Life Techonolgies, Paisley, Scotland), 60 nM carboxy-x-rhodamine (Cambridge Biosciences, Cambridge, UK), and 300nM of *UTR-R2* were prepared. One contained 300nM of *MBP-LR-ALL* and the other contained 300nM of *MBP-LR-1* (Oswel DNA Service, Southampton).

15 *MBP-LR-ALL* is a forward primer which is universal for all HCV genotypes and is derived from the primer *UTR-L2*. It has the following sequence:

5' FCA CCT TCA CCC TCA GAA GGM GCG UCT AGC CAT GGC GTT  
AG 3'

20

*MBP-LR-1* is a forward HCV type 1 specific primer which is derived from the primer *Spec1*. It has the following sequence:

5' FCA CCT TCA CCC TCA GAA GGM GCC GCT CAA TGC CTG GAG 3'

25

(F=FAM; M=MeREDdU and U=Uracil)

The results are shown in Figure 4 where the conditions are as follows:

	<b>1</b> (MBP-LR- 1 + UTR- R2)	<b>2</b> (MBP-LR- 1 + UTR- R2)	<b>3</b> (MBP-LR- 1 + UTR- R2)	<b>4</b> (MBP-LR- 1 + UTR- R2)	<b>5</b> (MBP-LR- 1 + UTR- R2)	<b>6</b> (MBP-LR- 1 + UTR- R2)	<b>7</b> (MBP- LR-1 + UTR- R2)	<b>8</b> (MBP-LR- 1 + UTR- R2)	<b>9</b> (MBP- LR-1 + UTR- R2)
<b>A</b>	1a	1a	1a	1a	3a	3a	-ve control	-ve control	-ve control